

CNME040405AP

Haemodialysis Powder A



The product is the special materials used for the preparation of haemodialysis dialysate whose function is removing metabolic waste and maintaining the balance of water, electrolyte and acid-base by the dialyser.

Description: white crystalline powder or granules

Application: The concentrate made from haemodialysis powder matching with haemodialysis machine is suitable for haemodialysis.

Specification: 1172.8g/person/bag; 2345.5g/2 person/bag; 11728g/10 person/bag

Dosage: 1 bag/patients; 1 bag/ 2 patients; 1 bag/ 10 patients

Storage: Sealed storage, avoiding direct sunlight, good ventilation and avoiding freezing should not be stored with the toxic, contaminated and bad smell goods.

Bacterial endotoxins: The product is diluted to dialysis by endotoxin testing water, bacteria endotoxins should not be more than 0.5EU/ml.

Insoluble particles: The product is diluted to dialysate, the particle content after deducting the solvent: 210um particles should not be more than 25's/ml; 225um particles should not be more than 3's/ml.

CNME040405BP

Hemodialysis Powder B



This product contains 84g sodium bicarbonate, and is the special materials used for the preparation of haemodialysis dialysate whose function is removing metabolic waste and maintaining the balance of water, electrolyte and acid-base by the dialyser.

Description: white crystalline powder or granules

Application: The concentrate made from haemodialysis powder matching with haemodialysis machine is suitable for haemodialysis.

Specification: 588g/person/bag; 1176g/2 person/bag; 5880g/10 person/bag

Dosage: 1 bag/patients; 1 bag/ 2 patients; 1 bag/ 10 patients

Storage: Sealed storage, avoiding direct sunlight, good ventilation and avoiding freezing, should not be stored with the toxic, contaminated and bad smell goods.

Bacterial endotoxins: The product is diluted to dialysis by endotoxin testing water, bacterial endotoxins should not be more than 0.5EU/ml.

Insoluble particles: The product is diluted to dialysate, the particle content after deducting the solvent: 210um particles should not be more than 25's/ml; 225um particles should not be more than 3's/ml.

CNME040405AF

Haemodialysis Fluid A



The product is the special materials used for the preparation of haemodialysis dialysate whose function is removing metabolic waste and maintaining the balance of water, electrolyte and acid-base by the dialyser.

Description: Dialysate A is colorless, clear and no precipitation liquid

Application: Haemodialysis Dialysate matching with haemodialysis machine is suitable for haemodialysis.

Specification: 5L/person/barrel; 10L/2 person/barrel; 11L/2 person/barrel;

Dosage: 5L of dialysate A for one patient; 10 of dialysate A for twopatient; 11L of dialysate A for two patient;

Storage: Sealed storage, avoiding direct sunlight, good ventilation and avoiding freezing, should not be stored with the toxic, contaminated and bad smell goods.

Bacterial endotoxins: The product is diluted to dialysis by endotoxin testing water, bacteria endotoxins should not be more than 0.5EU/ml.

Insoluble particles: The product is diluted to dialysate, the particle content after deducting the solvent: 210um particles should not be more than 25's/ml; 225um particles should not be more than 3's/ml.

CNME040405BF

Hemodialysis Fluid B



This product contains 84g sodium bicarbonate, and is the special materials used for the preparation of haemodialysis dialysate whose function is removing metabolic waste and maintaining the balance of water, electrolyte and acid-base by the dialyser.

Description: Dialysate B is colorless, clear and no precipitation liquid

Application: Haemodialysis Dialysate matching with haemodialysis machine is suitable for haemodialysis.

Specification: 7L/person/barrel; 10L/2 person/barrel; 12L/2 person/barrel;

Dosage: 7L of dialysate B for one patient; 10 of dialysate B for two patient; 11L of dialysate B for two patient;

Storage: Sealed storage, avoiding direct sunlight, good ventilation and avoiding freezing, should not be stored with the toxic, contaminated and bad smell goods.

Bacterial endotoxins: The product is diluted to dialysis by endotoxin testing water, bacterial endotoxins should not be more than 0.5EU/ml.

Insoluble particles: The product is diluted to dialysate, the particle content after deducting the solvent: 210um particles should not be more than 25's/ml; 225um particles should not be more than 3's/ml.

CNME040405BB

Concentrates For Haemodialysis (Dialysis Powder) Barrel



The product is the special materials used for the preparation of haemodialysis dialysate whose function is removing metabolic waste and maintaining the balance of water, electrolyte and acid-base by the dialyser.

Description: White crystalline powder or granules.

Application: The concentrate made from haemodialysis powder matching with haemodialysis machine is suitable for haemodialysis.

Specification: 650g/barrel; 700g/barrel; 720g/barrel.

Dosage: One barrel of product can be used to treat one patient.

Storage: Sealed storage, avoiding direct sunlight, good ventilation and avoiding freezing, should not be stored with the toxic, contaminated and bad smell goods.

Bacterial endotoxins: The product is diluted to dialysis by endotoxin testing water, bacteria endotoxins should not be more than 0.5 EU/ml.

Insoluble particles: The product is diluted to dialysate, the particle content after deducting the solvent: >10µm particles should not be more than 25's/ml; >25µm particles should not be more than 3's/ml.

Term of validity : 2 4 months from the date of production

CNME040405BB

Concentrates For Haemodialysis (Dialysis Powder) Poder Bag



The product is the special materials used for the preparation of haemodialysis dialysate whose function is removing metabolic waste and maintaining the balance of water, electrolyte and acid-base by the dialyser.

Description: White crystalline powder or granules.

Application: The concentrate made from haemodialysis powder matching with haemodialysis machine is suitable for haemodialysis.

Specification: 650g/bag; 700g/bag; 720g/bag

Dosage: One bag of product can be used to treat one patient.

Storage: Sealed storage, avoiding direct sunlight, good ventilation and avoiding freezing, should not be stored with the toxic, contaminated and bad smell goods.

Bacterial endotoxins: The product is diluted to dialysis by endotoxin testing water, bacteria endotoxins should not be more than 0.5 EU/ml.

Insoluble particles: The product is diluted to dialysate, the particle content after deducting the solvent: >10µm particles should not be more than 25's/ml; >25µm particles should not be more than 3's/ml.

Term of validity: 2 4 months from the date of production

CNME040405CA

Citric Acid DisinfectCart



Product Specification: 32g/Cartridge

Main Ingredient and Content: Citric Acid 100%

Scope: Combining thermal disinfection procedures removes calcium and magnesium deposits and sterilizes the water system inside the dialysis machine.

CNME040405CB

Sodium Carbonate CleanCart



Product Specification: 13g/Cartridge

Main Ingredient and Content: Sodium Carbonate 100%

Scope: Combining thermal disinfection procedures removes organic sediment, fat and protein of the internal water way system of dialysis machine.