



Product Certification

VMH-03-03-270 Issue Date 08-27-2024

Finished Good Code 5C4482		Batch No. H25H14138		Mfg. Date 2025-08-14		Sterilization Method: Gamma ☐ EO ☒		Sterilization Date 08-18-25	
Product Description MINICAP EXTENDED LIFE PD TRANSFER SET WITH TWIST CLAMP								Expiration 2030-08-14	
Biological Testing	Pyrogenicity Limit 0.5 EU/ml	☐ V13-01-K ☒ V130110313		Alternate Batch Testing ☐		Sample Result (Eu/ml): 0.0102		Date Passed: 08-27-25	
	Sterility (SAL 10 ⁻⁶)	Biological Indicators		Number tested: 16				Date 08-25-25	
Gamma Sterilization:									
Minimum KiloGray Dose per Spec.: NA					Maximum KiloGray Dose per Spec.: NA				
Minimum KiloGray Dose received.: NA					Maximum KiloGray Dose received.: NA				
☐ V20-03-03-027 ☐ V20-03-04-013 Final Physical Testing Ref: ☒ 20-03-20-012 ☐ V20-03-09-405 ☐ Other: _____								Results: PASS	
Comments and/or additional information (if needed)									
THIS LOT HAS BEEN PRODUCED IN ACCORDANCE WITH PROCEDURES IN THE MASTER PRODUCTION DOCUMENT.									
Plant of Manufacture: MOUNTAIN HOME, AR USA				Date Released 08-27-25		Quantity Released 12096 EA			
Prepared by/Date C. Demoupeul 8.29.25				QM Manager/Date Jon Abraham 9/2/25					

END OF FORM, ADDITIONAL PAGE(S) FOR INSTRUCTIONAL USE ONLY

Release Status: Issued - Not Yet Effective

INSTRUCTIONS TO COMPLETE FORM - VMHQCF2360

- A. Issue Date of SOP (VMH-03-03-270). Preparer will verify that current issue of form is being used.
- B. These are obtained from appropriate production router, BOM, AS400, POMS, and/or BRS.

- Code
- Batch
- Mfg Date (Manufacturing date is first day of production)
- Description
- Expiration Date
- Release Date
- Release Quantity

- C. This information is obtained from BRS or from VMH-QCF-1589 "Gas Sterilization and Radiation Sterilization Verification", if manual procedure was used.
- Sterilization Method
 - Sterilization Date (Date will be taken from the 1st load)

- D. Pyrogenicity -- Check box for the appropriate test method for batch.

- VMHQCF3988 (V130110313)
- V1190968 (V13-01-K)
- LIMS Printout

Sample results (Eu/ml) = Max Sample Conc. = Date passed

When LIMS printout is being used, the date entry is APPROVED/REVIEWED = Date Passed

NOTE: If pyrogen were tested using alternate testing, retrieve from "Alternate Batch Testing for Magellan or Apollo/Odyssey Pyrogen". If documented on two alternate batch testing dates, use both values and dates.

- E. Sterility – Enter number tested for the entire load and date passed. Retrieve information from "External Process Challenge Device Sampling and Testing Report" VMH-QCF-2018 and/or LIMS. Date Passed will be taken from last load processed.

When LIMS printout is used, the date entry is APPROVED/REVIEWED = Date Passed

Not applicable for Gamma sterilized product.

- F. Gamma Sterilization – Enter Min and Max Dose per Spec from Radiation Processing Parameter Sheet or PCS. Enter the Min and Max received from Certificate of Processing.

- G. Check box for appropriate Specification number from PCS or POMS. Verify conforming or passing results from Final Inspection Report, Quality Verification closure in POMS, or retention disposition.

- H. Individual that prepares form will sign/date.

- I. QM Manager or delegate approving certificate will sign/date.