

Haemodialysis

smartbag®

The smart way of concentrate supply



**FRESENIUS
MEDICAL CARE**

smartbag®

The environmentally-friendly flexible bag

The smartbag® – combining the advantages of a canister and a bag

Handling convenience

- Easy to carry
 - comfortable grip and reduced weight
- Stand-alone stability
- Compatible with all devices
 - equipped with standard suction rods

Environmentally friendly

- Free from PVC and plasticizers
- Recyclable packaging material¹
- 90% less disposal space needed
 - due to collapsible bag design²
- Up to 75% reduced plastic packaging material
 - in comparison to a conventional canister² due to a thin and lightweight bag



Compatible with all dialysis machines



Stand-alone stability



Reduced waste volume

The smartbag® – containing liquid acid concentrates (1+44 AC-F)

Compositions of ready to use dialysis fluid (after mixing with bicarbonate concentrate 8.4% and purified water)

Type	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ²⁺ mmol/L	Mg ²⁺ mmol/L	Cl ⁻ mmol/L	HCO ₃ ⁻ mmol/L	Acetate mmol/L	Glucose g/L	Osmolarity mosm/L	Art. No. 4.7 L Box of 60	Art. No. 4.7 L Box of 2
smartbag® 111.25	138.00	1.00	1.25	0,5	107.50	32.00	3.00	1.00	289	—	F00006883
smartbag® 111.50	138.00	1.00	1.50	0.5	108.00	32.00	3.00	1.00	290	F000 000 64	F000 006 4D
smartbag® 211.25	138.00	2.00	1.25	0.5	108.50	32.00	3.00	1.00	291	F000 000 65	F000 006 5D
smartbag® 211.50	138.00	2.00	1.50	0.5	109.00	32.00	3.00	1.00	292	F000 000 66	F000 006 6D
smartbag® 211.75	138.00	2.00	1.75	0.5	109.50	32.00	3.00	1.00	292	F000 002 32	F000 023 2D
smartbag® 311.25	138.00	3.00	1.25	0.5	109.50	32.00	3.00	1.00	293	F000 000 67	F000 006 7D
smartbag® 311.50	138.00	3.00	1.50	0.5	110.00	32.00	3.00	1.00	294	F000 000 68	F000 006 8D
smartbag® 311.75	138.00	3.00	1.75	0.5	110.50	32.00	3.00	1.00	295	—	F000 075 75
smartbag® 411.25	138.00	4.00	1.25	0.5	110.50	32.00	3.00	1.00	295	F000 002 36	F000 070 41
smartbag® 411.50	138.00	4.00	1.50	0.5	111.00	32.00	3.00	1.00	296	—	F000 075 76
smartbag® 411.75	138.00	4.00	1.75	0.5	111.50	32.00	3.00	1.00	297	—	F000 075 77

Acid concentrate per pallet

564 L

564 L

Bags/box – boxes/pallet

60–2

2–60

Support plate for smartbag (optional accessory)

for application with 5008 CorDiox

M442061

for application with 6008 CAREsystem, 5008S CorDiox

M423591

Code message: e.g. **smartbag® 111.50** (1xx.x = K⁺ in mmol/L) (x1x.x = Glucose in g/L) (xx1.5 = Ca²⁺ in mmol/L)

Technical changes reserved.

smartbag® 300L

The economic way to provide central concentrate supply

smartbag® 300L – facilitating central concentrate supply at reduced storage space

Handling convenience

- Space-saving in comparison to canisters
- Bag-in-Box concept making box opening obsolete
- Easy connection to the CDS3 Central Delivery System

Economic efficiency

- 600l of concentrate on one pallet
- Single-use packaging, eliminating the cost for return and cleaning of bulk containers
- No wastage of remaining concentrate

Environmentally-friendly product design

- Fully collapsible: Minimised waste volume
- 100% PVC-free
- Recyclable packaging material¹
- 32% improvement of the environmental performance in comparison to canisters³

Flexibility and convenience with the CDS3 Central Delivery System

- Modular design for adaptation to the building conditions of the dialysis unit
- Automatic switch to the standby smartbag® in case the smartbag® runs empty



Compatible with FME CDS3

The smartbag® – containing liquid acid concentrates (1+44 AC-F)

Compositions of ready to use dialysis fluid (after mixing with bicarbonate concentrate 8.4% and purified water)

Type	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ²⁺ mmol/L	Mg ²⁺ mmol/L	Cl ⁻ mmol/L	HCO ₃ ⁻ mmol/L	Acetate mmol/L	Glucose g/L	Osmolarity mosm/L	Art. No. 300 L
smartbag® 211.25	138.00	2.00	1.25	0.5	108.50	32.00	3.00	1.00	291	F000 011 68
smartbag® 211.50	138.00	2.00	1.50	0.5	109.00	32.00	3.00	1.00	292	F000 030 68
smartbag® 311.25	138.00	3.00	1.25	0.5	109.50	32.00	3.00	1.00	293	F000 030 69
smartbag® 311.50	138.00	3.00	1.50	0.5	110.00	32.00	3.00	1.00	294	F000 030 70
smartbag® 411.25	138.00	4.00	1.25	0.5	110.50	32.00	3.00	1.00	295	F000 066 59
Acid concentrate per pallet										600 L
Boxes / pallet										2

Code message: e. g. **smartbag® 211.50** (2xx.x = K⁺ in mmol/L) (x1x.x = Glucose in g/L) (xx1.5 = Ca²⁺ in mmol/L)

Technical changes reserved.



Bag-in-box system



Reduced waste volume



Easy connection and handling

References

- 1 According to Directive on Packaging and Packaging Waste (94/62/EC) and standard EN 13427:2004.
- 2 Comparison versus baseline product (rigid canister of 6 L 1+ 34).
- 3 Fresenius Medical Care (FME) internal study (2013): Comparative life cycle assessment of product group FME acid concentrates. Eco-performance is always calculated as average of 15 environmental impact categories versus the baseline product (rigid canister of 6 L 1+34).

smartbag is protected as trademark in selected countries

